

# CANDESCENCE® Foldover Blister

## A SUSTAINABLE OPTION FOR VISUAL PACKAGING



Designed for two-sided messaging and product visibility. Omni-directional adhesion for sealing and gluing without compromising your print surface. Our unique coating formulation delivers an ultra-smooth print surface on both the primary and secondary side for the most demanding applications.

### PRODUCT USES

- Visual packaging
- Computer cables and accessory packaging
- Personal electronics packaging
- Vitamin and health supplements packaging
- Nutraceutical packaging

### CANDESCENCE® FOLDOVER BLISTER TYPICAL VALUES

English	Units	14-176	16-195	18-212	20-230	24-263
Caliper (nominal)	mils	14.0	16.0	18.0	20.0	24.0
Basis Weight (nominal)	lbs/3msf	176	195	212	230	263
Brightness	GE	89	89	89	89	89
b* Value	Spectro	1.2	1.2	1.2	1.2	1.2
Coated Smoothness (primary side)	Sheffield	25	28	31	34	40
Parker Print (primary side)	microns	1.2	1.2	1.2	1.2	1.2
Plybond	lbs/in <sup>2</sup>	42	40	37	37	37
Stiffness – MD	g-cm	100	145	200	260	405
Stiffness – CD	g-cm	45	70	95	120	200
Metric	Units	356-286	406-317	457-345	508-374	610-428
Thickness	microns	356	406	457	508	610
Basis Weight	gsm	286	317	345	374	428
Stiffness – MD	mN	203	294	406	528	822
Stiffness – CD	mN	91	142	193	244	406